

# GluGo2.0

Assessing the Merits of the current GluGo Algorithm

# New T1D Paradigm on the Horizon



Dexcom G6 Sensor



**NIGHTSCOUT**  
#WeAreNotWaiting

# GluGo's Algorithm

- Structural model (no machine learning involved)
- Jingxian Liu (17') and Jihong Jin (16')

$$UpwardSlope(/min) = \frac{\frac{Sens}{CarbRatio} * Carb}{4 * 60}$$

$$DownwardSlope(/min) = -\frac{Insulin * Sens - Basal * 4 * Sens}{4 * 60}$$

# Initial Hypotheses, Questions, and Concerns

- **Hypotheses**

- Algorithm will be fairly accurate overall
- substantially less accurate at predicting high glucose

# Initial Hypotheses, Questions, and Concerns

- **Hypotheses**

- Algorithm will be fairly accurate overall
- substantially less accurate at predicting high glucose

- **Questions**

- Are the predictions biased?
- Why is it less accurate for high glucose


# Initial Hypotheses, Questions, and Concerns

- **Concerns**

- The algorithm is based solely on theory, and makes ***many*** simplifying assumptions
- Doesn't "learn" about the patient (no ML) – not personalized

1h 1d 1w 1m YTD 1y all

# Predicted vs. Actual Glucose Over 1 24 Hour



# Summary Statistics on Residuals and Absolute Residuals

## Residuals:

Mean	Median	Standard Deviation	Minimum	Maximum
-2.92	-2	9.11	-285	296


Sample Size
n = 28,078

## Absolute Residuals:


Mean	Median	Standard Deviation	Minimum	Maximum	Average % Residual
5.16	3	8.05	0	296	4.01%



**Residuals (Predicted – Actual) by  
Actual Glucose Value**



**Residual Density Plot**





# Upon Closer Inspection...

- The algorithm predicts at most **15 minutes** into the future
  - Usually only **5 minutes**
- The test data used by the previous group is **not a random sample**
  - Only a fraction of the total sample includes carb intake
- Algorithm functions in two “modes”
  - **Mode 1**: After Carbohydrate Intake
  - **Mode 2**: No Carbohydrate Intake

# Mode 1 (Carbs) vs. Mode 2 (No Carbs)

	Mode 1	Mode 2				
Mean Residual	-5.53	-3.25				
Mean Abs. Residual	9.74	4.23				
Average % Difference	7.37%	3.73%				
Sample Size	n = 3,016	n = 1,714				

# Mode 1 (Carbs) vs. Mode 2 (No Carbs)



# Conclusions about Algorithm

1. Old algorithm in it's current form is not optimal

1. Prior Assessments of the model's accuracy are inflated due to a biased sampling procedure

# Looking Forward: Building New Test Datasets

- Compiling a new dataset

A	B	C	D	E	F	G
Timestamp	Glucose (mg/dL)	Basal Insulin (U/hr)	Insulin Sensitivity (mg/dL/U)	Carb Ratio (g/U)	Bolus Insulin (U)	Carb Intake (g)

# Looking Forward: Building New Test Datasets

- Compiling a new dataset

A	B	C	D	E	F	G
Timestamp	Glucose (mg/dL)	Basal Insulin (U/hr)	Insulin Sensitivity (mg/dL/U)	Carb Ratio (g/U)	Bolus Insulin (U)	Carb Intake (g)

- Testing the new dataset

```

from dataAccess import dataAccess as DB
import datetime as DT
import csv

#path to filled out DB
dbpath = "lab.db"
##Connect to the Database by creating an instance of dataAccess
katDB = DB("Kate",dbpath)

for patient in [katDB]:
    dat = patient.getData()

gluPredict = []

i = 0
while i < len(dat):
    # if dat[i][3] == None:
    #     gluPredict.append()
    if dat[i][9] == "INRANGE":
        if dat[i][3] != None:
            bolis = dat[i][3]
            carb = dat[i][4]
            carbRate = dat[i][15]
            basalRate = dat[i][16]
            sens = dat[i][17]
            gluPredict.append(dat[i-1][2])


            currentTime = dat[i][1]
            # Evan Laferriere (a year ago) · Add files via upload
            upwardSlope = ((sens/carbRate) * carb) / (4*60.0)
            downwardSlope = (- bolis * sens - basalRate * 4 * sens) / (4*60.0)
            i = i+1

        if dat[i][3] == None:
            upValue1 = round(dat[i-1][2] + ((dat[i][1] - dat[i-1][1]).total_seconds() / 60.0) * upwardSlope)
            upValue2 = round(upValue1 + ((dat[i+1][1] - dat[i][1]).total_seconds() / 60.0) * upwardSlope)
            upValue3 = round(upValue2 + ((dat[i+2][1] - dat[i+1][1]).total_seconds() / 60.0) * upwardSlope)
            i = i+3
            gluPredict.append(upValue1)
            gluPredict.append(upValue2)
            gluPredict.append(upValue3)

        while dat[i][3] == None and i < len(dat) -4 and dat[i+1][9] == "INRANGE" and (dat[i][1] - currentTime).total_seconds() < 15300 :
            newValue = round(dat[i-1][2] + ((dat[i][1] - dat[i-1][1]).total_seconds() / 60.0) * downwardSlope)
            gluPredict.append(newValue)
            i = i+1

```

# Moving from a database to directly reading a csv and using pandas library



```

import csv
import pandas
import numpy

data = pandas.read_csv('mlk.csv')

#for i in range(0,3):
#     print(data["Timestamp"][i])    Testing if pandas work.

length_csv = len(data)
#print(length_csv)
gluPredict = []

i = 0
while i < length_csv:
    #print(data["Bolus Insulin (U)"][i]) just to check the values.
    # if dat[i][3] == None: LEGACY CODE
    #     gluPredict.append() LEGACY CODE
    if data["Bolus Insulin (U)"][i] != None : # assuming that if Bolus insulin has a value then, it will be an
        bolis = data["Bolus Insulin (U)"][i]
        carb = data["Carb Intake (g)"][i]
        carbRate = data["Carb Ratio (g/U)"][i]
        basalRate = data["Basal Insulin (mg/dL)"][i]
        sens = data["Insulin Sensitivity (mg/dL/U)"][i]
        gluPredict.append(data["Glucose (mL/dL)"][i])

        print("glucose level: ")
        print(data["Glucose (mL/dL)"][i]) # print the glucose level
        currentTime = data["Timestamp"][i] # put it as a time stamp but its not really time stamp. Not sure what the value means in their database.

        upwardSlope = ((sens/carbRate) * carb) / (4*60.0)
        downwardSlope = (- bolis * sens - basalRate * 4 * sens) / (4*60.0)
        i = i+1
        if data["Bolus Insulin (U)"][i].isnull():
            upValue1 = round(data["Glucose (mL/dL)"][i-1] + ((data["Timestamp"][i] - data["Timestamp"][i-1]).total_seconds() / 60.0) * upwardSlope)
            upValue2 = round(upValue1 + ((data["Timestamp"][i+1] - data["Timestamp"][i]).total_seconds() / 60.0) * upwardSlope)
            upValue3 = round(upValue2 + ((data["Timestamp"][i+2] - data["Timestamp"][i+1]).total_seconds() / 60.0) * upwardSlope)
            i = i+3
            gluPredict.append(upValue1)
            gluPredict.append(upValue2)
            gluPredict.append(upValue3)

        while data["Bolus Insulin (U)"][i] == None and i < length_csv -4 and (data["Timestamp"][i] - currentTime).total_seconds() < 15300 :
            newValue = round (data["Glucose (mL/dL)"][i-1] + ((data["Timestamp"][i] - data["Timestamp"][i-1]).total_seconds() / 60.0) * downwardSlope)
            print("i am here")
            gluPredict.append(newValue)
            i = i+1

    else:
        newValue = round(data["Glucose (mL/dL)"][i-1] - ((data["Timestamp"][i] - data["Timestamp"][i-1]).total_seconds() / 60.0) * (data["Basal Insulin"][i] *

```



# Looking Forward: Better Algorithms

- SVR Model with Physiological Features:

